

## CM50170-01 Oceanographic CO<sub>2</sub> Coulometer Dissolved Inorganic Carbon (DIC) in Sea Water For use with CM5260 MODICA system



The Coulometer uses coulometric detection. The carbon Coulometer measures carbon as CO<sub>2</sub>. The gas stream resulting from the MODICA unit is bubbled through the coulometer analytical cell. The carbon coulometer solution contains ethanolamine and a colorimetric pH indicator. The CO<sub>2</sub> from the gas stream reacts with the ethanolamine forming a strong titratable acid, causing the color indicator to fade. The coulometer photometer recognizes this condition and initiates the electrochemical generation of electrons returning the solution to the original color. The current for this 100% efficient coulometric process is integrated and digitally displayed in user selected units. The Oceanographic unit is supplied with the CM5011 Emulation firmware necessary to allow direct connection to the MODICA software. No other software is required.


The CM50170 Coulometer sits on a conventional lab bench capable of supporting 40 lbs. The Coulometer stands 12" wide by 12" high by 19.2" deep. The instrument is constructed of aluminum and steel. The unit is designed with a cell compartment, power switch, cell current switch, and a 13.3 x 12.1" LCD touch screen to act as the user interface to the instrument. The instrument is supplied with is an analytical cell assembly, power conditioner, power cord with a NEMA 5-15 plug and an RS232 serial cable. The unit is also supplied with a set of cell reagents and operation manual. The unit when supplied for 110-120V 50/60 HZ operation requires one AC circuit capable of supplying 1.5 amps.

### CM50170 CO<sub>2</sub> Coulometer

- No user calibration
- 100% efficient coulometric detection
- Wide, linear dynamic range
- Readability to 0.01 µg Carbon
- Relative standard deviations of < 0.2% for standard certified materials
- 12.1" LCD Touch Screen
- Typical analysis time of 8-10 minutes
- LIMS Compatible

#### Part Numbers

CM5017-01 for 110V / 50/60Hz  
CM5017-02 for 220V / 50/60Hz




## CM210-032 – Titration Cell includes:

1. CM200-051 – Titration Cell with Side Arm
2. CM119-027 – Cathode Top, White Teflon
3. CM101-135 – Platinum Electrode (black lead)
4. CM190-020 – Gas Inlet Tube (blue tag)
5. CM190-021 – Gas Exit Tube (red tag)
6. CM119-028 – Anode Top, White Teflon
7. CM101-033 – Silver Electrode (red lead)
8. CM121-001 – Stir Bar

## CM310-001 – Cell Reagent Kit includes:

- CM300-001 – Carbon Cathode Solution (1 gallon)
- CM300-002 – Carbon Anode Solution (16 ounces)
- CM300-003 – Potassium Iodide (50 grams)



## Technical Differences between CM50170 and CM50150 Oceanographic CO<sub>2</sub> Coulometers

1. Touch screen - New touch screen streamlined for Oceanographic application.
2. Emulation Standard - CM50170 is supplied with CM5011 Emulation as standard and only program.
3. 50 mA Optimization – Instrument is supplied with 50 mA cell current. Instrument is optimized so that the full analytical measuring range of the A to D converter is used for the 50mA cell current. Instrument can only be used with a 50-mA cell current.
4. Open Cell Compartment- Instrument is supplied with an external open style cell compartment and associated Cell assembly for ease of use and compatibility to older style CM5011 configuration.